
PRODIGY

AutoCode[®]

Blood Glucose Monitoring System

Owner's Manual

Speaks in
4 Languages

PRODIGY AutoCode® Blood Glucose Meter

Dear Prodigy AutoCode® Meter Owner:

Thank you for choosing the Prodigy AutoCode® as your blood glucose monitoring system. Please read this manual carefully as it contains important information about your new Prodigy AutoCode® system.

The Prodigy AutoCode® meter has the latest technology for blood glucose monitoring and is made with quality components. It is easy to use, provides fast test results with a minimal sample of blood, has a large, easy-to-read display screen, and is lightweight and portable for your convenience.

The Prodigy AutoCode® meter is designed to help you and your healthcare professionals monitor your blood glucose levels. This Owner's Manual will help you learn how to use your Prodigy AutoCode® meter effectively.

A warranty registration card is included with your system. Please return the completed card to us, complete the warranty online at www.prodigymeter.com, or call Customer Care to register your Prodigy AutoCode® meter. Customer Care is available by calling 1.800.243.2636 Monday - Friday 8:30AM to 5:30PM (EST). All questions concerning test results should be directed to your healthcare professional.

Prodigy AutoCode® meters:

- Require No Coding, allowing you to save time and decreases the chances of human error due to improper coding.
- Allow you to perform Alternate Site Testing (AST).
- Have memory and data management capabilities.

Prodigy's Diabetes Management software gives you and your healthcare professional, powerful graphic tools to manage your diabetes.

Table of Contents

Letter to Meter Owner	2-3
Table of Contents	4-5
Important Safety Instructions	6-9
Important Health-Related Information	10
Alternate Site Testing (AST)	11-12
Introduction to Prodigy AutoCode® Meters	13-14
Intended Use	13
Test Principle	13-14
About Prodigy No Coding Blood Glucose Test Strips	14-16
Key Functions of the Test Strip	14-15
Important Test Strip Information	15-16
Contents of the Prodigy AutoCode® Meter	17
Setting up the Prodigy AutoCode® Meter	18-27
Key Functions	18
Meter Display	19-20
Speaking Function	21
Setup Steps	22-27
Control Solution Testing	28-34
Performing a Control Solution Test	31-34
Performing a Blood Glucose Test	35-44
Preparing the Lancing Device	35-36
Getting a Blood Sample	37-38
Tips for Applying Blood Without Sight	39-41
Applying Blood to the Test Strip	42-44

Table of Contents

Using the Meter Memory	45-47
Viewing Results on a Personal Computer	47
Specifications	48
Caring for your Prodigy AutoCode® Meter	49-54
Cleaning your System	50-51
Lancing Device.....	51
Disinfecting your System.....	51-52
Storing your System	52
Battery	53
Low Battery.....	53
Battery Replacement.....	54
System Troubleshooting	55-61
Special Messages	55-56
Error Messages.....	56-57
Troubleshooting Guide	58-61
Information About your Prodigy AutoCode® Meter	62-71
Comparing Meter and Laboratory Results	62-63
Performance Characteristics.....	64-65
Expected Test Results.....	66
Unexpected Test Results	67-68
Accuracy of Meter Measurements	69-71
Symbol Information	72
Warranty Information	73
Limited Lifetime Warranty	73

Important Safety Instructions

READ ALL INSTRUCTIONS THOROUGHLY BEFORE USING.

The following basic safety precautions should always be taken.

- Close supervision is necessary when the device is used by, on, or near children and handicapped persons.
- Use the device only as described in this manual.
- Use only Prodigy No Coding Blood Glucose Test Strips and Prodigy Control Solution with your Prodigy AutoCode® meter. Using other test strips and control solutions with this meter can produce inaccurate results.
- Do not use the device if it has suffered any damages or if it appears to not be working properly.
- Keep the test strip vial away from children; the vial cap and test strips can be a potential choking hazard.
- Never chew or swallow a test strip. If this occurs, seek medical assistance immediately.
- All parts of the system are considered biohazardous and can potentially transmit infectious diseases, even after cleaning and disinfection.

Important Safety Instructions

- The device is intended for patient over-the-counter use to self-monitor blood glucose. It should not be used to collect blood from more than one person.
- The meter and lancing device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients!
- To reduce the risk of electromagnetic interference, **DO NOT** use this meter near cellular and cordless telephones, intercom systems, electric garage gates, radio transmitters and other electrical or electrical equipment that are sources of electromagnetic radiation, which may interfere with the meter's operation.

For further information, please refer to:

FDA Public Health Notification: Use of Fingerstick Devices on More Than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication (2010) <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens (2011)
<http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>

DO NOT CHANGE YOUR TREATMENT BASED ON A SINGLE RESULT THAT DOES NOT MATCH WITH HOW YOU FEEL, OR IF YOU BELIEVE THAT YOUR TEST RESULT COULD BE INACCURATE.

Limitations

Hematocrit: Hematocrit levels below 20% or above 60% can cause false results. Please consult your healthcare professional if you do not know your hematocrit level.

For self-testing only.

For in vitro diagnostic use only.

Do not share your Prodigy AutoCode® meter, test strips, lancets and lancing device with anyone, including other family members.

Do not use on multiple patients.

The system should NOT be used for the diagnosis of or screening for diabetes.

The system should NOT be used in neonates.

The system should NOT be used on critically ill patients.

Not for use on patients who are dehydrated, hypotensive, in shock or in a hyperosmolar state.

Alternate Site Testing results should NOT be used for calibrating CGMs or for insulin dosing calculations when glucose levels may be rapidly changing or fluctuating, such as after meals or exercise, or after taking insulin.

Do not use the Prodigy AutoCode® Blood Glucose Monitoring System during or after a Xylose Absorption Test. This may falsely raise glucose results. Please check with your Doctor before using the Prodigy AutoCode® Blood Glucose Monitoring System.

The system exhibits interference from Acetaminophen, ascorbic acid and Ibuprofen. If you are taking acetaminophen containing drugs (Tylenol, etc., >8.0 mg/dL in your blood), or ibuprofen containing drugs (Advil, etc., > 50.0 mg/dL in your blood) or vitamin C (ascorbic acid, > 5.0 mg/dL in your blood), you may obtain inaccurate readings from this blood glucose monitoring system. Please check with your doctor before using the Prodigy Autocode® Blood Glucose Monitoring System.

High uric acid concentrations (>8 mg/dL in your blood) may interfere with the glucose measurements with this device. If you have medical conditions that are associated with high uric acid level or hyperuricemia (e.g. gout), then please check with your Doctor before using the Prodigy AutoCode® Blood Glucose Monitoring System.

Important Health-Related Information

Severe dehydration and excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult with a healthcare professional immediately.

If you are experiencing symptoms that are not consistent with your blood glucose test results and you have followed all instructions described in this owner's manual, contact your healthcare professional.

Test results below 60 mg/dL (3.3 mmol/L) indicate low blood glucose (hypoglycemia).

Test results greater than 240 mg/dL (13.3 mmol/L) indicate high blood glucose (hyperglycemia).

If your test result is below 60 mg/dL or above 240 mg/dL, repeat the test, and if the result is still below 60 mg/dL (3.3 mmol/L) or above 240 mg/dL (13.3 mmol/L), consult with your healthcare professional immediately.

Inaccurate results may occur in severely hypotensive individuals, individuals in shock, and in those experiencing a hyperglycemic-hyperosmolar state. Please refer to the test strip package insert for additional important information.

Alternate Site Testing

There are important limitations to performing AST. Consult with your healthcare professional before testing with a blood sample obtained from AST locations.

What is AST?

Alternate Site Testing (AST) means you can use parts of the body other than your fingertips to check your blood glucose levels. The system allows you to test from the palm, forearm, upper arm, calf or thigh, with equivalent results to fingertip testing.

What is the advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). At other body sites, nerve endings are not so numerous and you will not feel as much pain as you will experience at the fingertip.


Figure 1

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at fingertip reflects these changes faster than capillary blood at other sites. Therefore, if you are testing blood glucose levels during or immediately after a meal, physical exercise or stressful event, take the blood sample from your fingertip only.

Use AST only:

- In a pre-meal or fasting state (more than 2 hours after the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.
- During steady state blood glucose conditions.

Do NOT use AST if:

1. You have reason to believe you have hypoglycemia or hyperglycemia.
2. Your routine glucose results are often fluctuating.

Introduction to Prodigy AutoCode® Meters

Intended Use

The Prodigy AutoCode® Blood Glucose Monitoring System is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips, forearm, upper arm, palm, calf or thigh.

The system is intended for self-testing outside the body (in vitro diagnostic use) by people with diabetes and is an at-home aid to monitor the effectiveness of diabetes control. It should not be used for the diagnosis or screening of diabetes or for neonatal use.

Test Principle

Blood glucose is measured by an electrical current that is produced when a blood sample mixes with the reagent (special chemicals) of the test strip. The electrical current changes with the amount of glucose in the blood sample. The Prodigy AutoCode® meter measures the strength of the electrical current, calculates the blood glucose level and then displays the result in milligrams per deciliter (mg/dL) or millimoles per liter (mmol/L). (USA standard: mg/dL)

PRODIGY AutoCode® Blood Glucose Meter

The Prodigy AutoCode® meter, test strips and control solution has been designed, tested and is proven to work together as a system to produce accurate blood glucose test results.


Important:

Use only Prodigy No Coding Blood Glucose Test Strips and Prodigy Control Solution with your Prodigy AutoCode® meter. Using other test strips and control solutions with this meter can produce inaccurate results.

About Prodigy No Coding Blood Glucose Test Strips

Key Functions of the Test Strip

The Prodigy AutoCode® meter measures the amount of sugar in the blood, commonly referred to as blood glucose. Blood is applied to the opening of the absorbent channel on the end of the test strip and is automatically drawn into the test strip.

- 1. Absorbent Channel:** Apply a drop of blood here.
- 2. Confirmation Window:** To confirm if enough blood is applied or has been drawn into the strip.
- 3. Test Strip Handle:** Hold this part to insert the test strip into the test strip slot.
- 4. Contact Bar:** Insert this end into the test strip slot.

Please refer to the "Performing a Blood Glucose Test" section for complete instructions.

Important Test Strip Information

- Store test strip packages in a cool, dry place between 39.2°F-104°F (4°C-40°C) and 10-85% relative humidity.
- Do not freeze.
- Keep away from direct sunlight and heat. (Do not store in your car.)
- Store your test strips in their original vial only. Do not transfer them to a new vial or any other container.

PRODIGY AutoCode® Blood Glucose Meter

- Close the test strip's vial lid immediately after removing a strip, and use each strip promptly.
- With clean, dry hands you may touch the test strip anywhere on its surface when removing it from the vial or inserting it into the meter.
- Wash the test site thoroughly with soap and water.
- Record the discard date on the vial label when you first open it. Discard remaining test strips 90 days after the first opening date.
- Do not use test strips beyond the expiration date printed on the package.
- After removing the test strip from the vial, insert it directly into the meter. Do not lay the test strip on any other surface.

Warning: Keep the test strip vial away from children; the vial cap and the test strips can be a potential choking hazard. Never chew or swallow a test strip. If this occurs, please seek medical assistance immediately.

Contents of the Prodigy AutoCode® Meter

Please review the contents of your purchase to confirm that all the components are included as listed on the side of your meter box.

You will need the following items to test your blood glucose levels:

1. Prodigy AutoCode® Blood Glucose Meter
2. Prodigy No Coding Blood Glucose Test Strips
3. Lancets
4. Lancing Device
5. Prodigy Control Solution

If any item is not included in the box, please purchase it separately. Please contact Customer Care at 1.800.243.2636.

Setting Up the Prodigy AutoCode® Meter

Key Functions

Test Strip Slot

Insert the test strip here.
The meter will turn on automatically.

LCD Display Screen

Shows test results and messages.

M Button

Turns the meter on and performs other functions described in this manual.

USB Data Port

USB cable connection to your computer.

Set Button

Used to set up the meter. Located on the back of the meter, inside the battery compartment.



Figure 4

Meter Display

Meter Display

Day Average

Indicates that the displayed test result is an average.


Figure 5

Memory Symbol

Appears when reviewing the memory.

Temperature Symbol

Appears when temperature is outside of the operating range.

Low Battery Symbol

Appears when the battery power is low.

Voice Symbol

Shows if audio function is on.

Speaking Function

The Prodigy AutoCode® meter talks you through the testing procedure.

When does the meter speak?	What does the meter say?
When the meter is turned on.	"Your meter is on. Please wait. Please insert the test strip."
When the room temperature is outside of the operating range; 50°F-104°F (10°C-40°C).	"Temperature is out of range."
When the meter is ready to test (💧 appears on display).	"Please apply blood into the test strip."
When the test is completed (the result appears on the display).	"Blood glucose (number) milligrams per deciliter/ millimoles per liter."
When you turn off the meter.	"Goodbye."
When a used test strip is inserted.	"Test strip has been used."
When blood is applied to the test strip.	"Testing."

Setup Steps

The Prodigy AutoCode® meter has several different settings that may need to be adjusted before use. You should check the time and date each time the batteries are replaced.

STEP 1: Insert the Batteries

Open the battery cover located on the back of the meter. Insert two (2) AAA Alkaline batteries and align them properly. The meter will turn on and enter Setting Mode.

STEP 2: Enter Setting Mode

If your meter is off, press the **"SET"** button located in the battery compartment. The meter will enter setting mode.

Voice Speaks: **"Your meter is on.
Please wait. Setting mode. English."**

STEP 3: Select Language

The language option appears first with **"L-1"** on the LCD. **"L-1"** indicates the default language (English); **"L-2"** indicates the second language (Spanish); **"L-3"** indicates the third language (French); and **"L-4"** indicates the fourth language (Arabic). See figure 6.

Press the "M" button to select a language. With the correct language selected, press the "SET" button and a number will appear.

Figure 6

STEP 4: Set the Volume Level

- Numbers 1 to 3 indicate the speakervolume from low to high. If the speaker function is set on one of these levels, "🔊" will be present on the display during testing. See Figure 7.
- Number 0 indicates that the speaking function is turned off. See Figure 8.

Figure 7

Figure 8

Press the "M" button to select the speaking volume, and then press the "SET" button. The year segment will appear and start flashing. See Figure 9.


Figure 9

STEP 5: Set the Year

Press and release the "M" button to advance the year. With the correct year on the display, press the "SET" button. The date will then appear on the display with the month segment flashing. See Figure 10.


Figure 10

STEP 6: Set the Month

Press and release the "M" button to advance the month. With the correct month on the display, press the "SET" button. The date segment will start flashing. See Figure 11.


Figure 11

STEP 7: Set the Date

Press and release the "M" button to advance the date. With the correct date on the display, press the "SET" button. The time will appear on the display with the hour segment flashing. See Figure 12.


Figure 12

STEP 8: Set the Hour

Press and release the "M" button to advance the hour. With the correct hour on the display, press the "SET" button. The minutes segment will start flashing. See Figure 13.


Figure 13

STEP 9: Set the Minutes

Press and release the "M" button to advance the minutes. With the correct minutes on the display, press the "SET" button and then the current unit of measurement will start flashing.

Important: Day averages are calculated from results obtained during the 7, 14 and 28 days preceding the current date and time settings.

STEP 10: Select Unit of Measurement

Press and hold the “**M**” button for four (4) seconds until the unit of measurement you are choosing appears on the display. See Figure 14 and Figure 15. To skip this step, press the “**SET**” button.

Your meter can display test results in milligrams per deciliter (mg/dL) or millimoles per liter (mmol/L).

- The mg/dL unit is standard in the United States.
- The mmol/L unit is standard in Canada.
- The Prodigy AutoCode® meter is set by default to read glucose values in mg/dL, unless changed by the user.
- Reading glucose values in the incorrect units may cause serious problems such as improper dose of medication.


Figure 14


Figure 15

STEP 11: Delete Memory

When the “**dEL**” and flashing “**M**” symbols appear on the display, you can choose to clear the memory. If you do not want to clear the memory, press the “**SET**” button to skip this step.

If you want to clear ALL memory, press the “**M**” button for four (4) seconds. See Figure 16.

“**OK**” will appear to show that all memory has been deleted.


Figure 16

Setup is now complete. “**OFF**” will display and the meter will turn off. See Figure 17.

Figure 17

Control Solution Testing

Prodigy Control Solution should come with an insert.

IVD

• For in vitro diagnostic (external) use only.

• Please review all instructions before using.

INTENDED USE:

Prodigy Control Solution is used to check the testing technique and system performance. The Prodigy AutoCode® Blood Glucose Monitoring System can use both Prodigy Control Solution Low (Level 1) and High (Level 2).

SUMMARY:

Prodigy Control Solution is a red liquid that contains glucose. The control solution reacts with the test strips to produce a test result. When the control solution test result falls within the control range given for the lot of test strips being used, this is indication that the system is working as it should and your testing technique is good. The control range can be found on the test strip vial label.

CAUTION: The control solution ranges printed on the back of the Prodigy test strip vials are not blood glucose target ranges.

IMPORTANT INFORMATION:

- Prodigy Control Solution is only to be used with Prodigy Blood Glucose Monitoring Systems.
- Check the expiration date on the bottle. Do not use if expired.
- Use within 90 days after opening. Record the discard date on the bottle to remind you when to discard your control solution.

WHY PERFORM A CONTROL SOLUTION TEST:

- To ensure that your meter and test strips are working properly together.
- To allow you to practice testing without using your own blood.

WHEN TO PERFORM A CONTROL SOLUTION TEST:

- When you begin using a new vial of test strips.
- At least once a week.
- When you suspect the meter or test strips are not working properly.
- When your blood glucose test results are not consistent with how you feel, or when you think your test results are inaccurate.

PRODIGY AutoCode® Blood Glucose Meter

- If you drop the meter.
- If you change the batteries.

STORAGE AND HANDLING:

- Storage temperatures range: 35.6°F to 86°F (2°C to 30°C).
- Do not freeze.

COMPOSITION:

INGREDIENT*	LOW (Level 1)	HIGH (Level 2)
Glucose	0.05	0.19
Other non-active ingredients	99.95	99.81

*w/v %

TEST PROCEDURE:

- Check the expiration and discard dates on your control solution and test strip vial.
- Check your control solution to see if it is a Low (Level 1) or High (Level 2) solution.
- Control solution should come to room temperature (68°F to 77°F / 20°C to 25°C) before testing.

Performing a Control Solution Test

Start with the meter off.

STEP 1: Wash your Hands

Wash your hands with mild soap and water before performing any test. Be sure to dry them thoroughly.

STEP 2: Insert Test Strip

Remove a new test strip from the test strip vial, making sure to tightly replace the vial cap. Immediately insert the test strip into the meter's test strip slot. The meter will turn on automatically. When the "  " symbol is blinking, it means the meter is ready for testing.

STEP 3: Mark as a Control Solution Test

After the "  " symbol appears on the display, press the "M" button and a "  " symbol will appear on the display. With the "  " symbol on the display, the meter will not store your test result in the memory. If you decide not to perform a control solution test, press the "M" button again and the "  " symbol will disappear.

Important: Make sure that the meter is in control solution mode so the test result will not be stored in the meter's memory.

STEP 4: Apply Control Solution

1. Shake the control solution bottle well, then remove the cap.
2. Squeeze the bottle and discard the first drop, then wipe the dispenser tip with clean tissue paper or cotton.
3. Squeeze the bottle again to get a second drop onto a clean, non-absorbent surface or on your finger. See Figure 18.
4. Bring the tip of the test strip to touch the drop of solution until the meter beeps.


Figure 18

Caution: To avoid contaminating the control solution with the content of the test strip, do not apply the control solution directly onto the test strip from the bottle. See Figures 19 and 20.


Figure 19


Figure 20

STEP 5: Check if the Test Result is in Range.

After the meter counts down from 6 to 1, the control solution test result will appear on the display. Compare this test result with the range (Low-Level 1 or High-Level 2) on the back of the test strip vial. The control range can change with each new vial of test strips. Always use the control solution range on the label of your current vial of Prodigy No Coding Blood Glucose Test Strips. The result should fall within the printed range on the test strip vial. See Figure 21.


Figure 21

QUESTIONABLE RESULTS:

If the test result falls outside the range printed on the test strip vial, check the “Troubleshooting Guide” located in the “System Troubleshooting” section of the Owner’s Manual and repeat the test.

Results that fall outside of the control range may indicate:

- An error in how the test was performed
- Control solution bottle not shaken well
- Expired or contaminated control solution
- Control solution that is too warm or too cold
- Test strip deterioration
- Meter malfunction
- Test strips expired

Caution: Do not use the test strips or control solution after the expiration date, as the results may be inaccurate.

Contact Customer Care

If your test strips have been damaged or if you continue to get out-of-range results, it means that the system or control solution may not be working properly.

Do NOT use the system to test your blood glucose level and contact Customer Care at 1.800.243.2636.

Performing a Blood Glucose Test

Warning:

1. Lancets are for single use only.
2. Keep lancing device and lancets clean.
3. Use caution when removing the used lancet from the device and when disposing of it.
4. The meter and lancing device are for single patient use. Do NOT share them with anyone including other family members. Do NOT use on multiple patients.

Preparing the Lancing Device

STEP 1: Twist off the Lancing Device Cover. See Figure 22.


Figure 22

STEP 2: Insert a sterile lancet into the Lancet Carrier. See Figure 23.


Figure 23

STEP 3: Hold the Lancet Needle Cover and gently twist it off the Lancet. See Figure 24.


Figure 24

STEP 4: Re-cap the Lancing Device Cover.

STEP 5: The Adjustable Tip offers 10 levels of skin penetration. To select the desired depth, twist the Adjustable Tip in either direction until the preferred setting lines up with the arrow. See Figure 25.

To select the best depth:

1 - 3 for soft or thin skin

4 - 6 for average skin

7 - 10 for thick or calloused skin



Figure 25

STEP 6: Slowly pull the Loading Barrel away from the Lancing Device Cover until you hear a click, indicating that the Lancet Carrier is locked into position. Release the Loading Barrel to return it to its original position. See Figure 26.


Figure 26

The Lancing Device is prepared and ready to lance for a blood sample.


Getting a Blood Sample

STEP 1: Wash your hands and the puncture site with warm, soapy water. Rinse and dry completely. Warm your fingers to increase blood flow.

STEP 2: Select and Lance a Puncture Site.

- **Fingertip**

Position the end of the Adjustable Tip against the side-tip of the finger or alternate testing site. Press the Trigger Button, and then lift the Lancing Device away from the testing site after the puncture is complete.

- **For AST**


Please refer to the "Alternate Site Testing (AST)" section. Consult with your healthcare professional before performing AST.

STEP 3: Gently massage your finger or the puncture site to obtain the required blood volume. To perform the test, you need 0.7 μL of blood. For best results, wipe off the first drop of blood and gently massage the area to obtain another drop of blood. Do not smear the blood sample.

Step 4: Remove the Lancing Device Cover and push the Lancet Ejector to release the Lancet. Dispose of the used Lancet into the appropriate Sharps or Biohazard container.

Note: The first drop of blood usually contains tissue fluid and serum which may affect the test result. It should be discarded.

Tips for Applying Blood Without Sight

If you have little or no vision, using the tips below may help you successfully place a blood sample into the test strip.

Other techniques may work as well, but the following examples illustrate ways to bring the sample and the test strip together. Keep experimenting until you find a technique that consistently works for you.

Tip One - Smaller has its Advantages

Because the Prodigy AutoCode® meter is easy to move around, you can experiment with many different approaches to determine which works best for you. Get comfortable holding the meter in one hand while testing. The meter can be moved towards your finger as easily as your finger can be moved to the test strip.

Tip Two - Aim for the Tip

The tip of the test strip contains an opening to the absorbent channel where the blood is applied. When the tip touches the blood sample, the absorbent channel draws the blood into the test strip, much like a straw. Blood should be applied only on the very tip. Blood that is applied to the bottom or sides of the test strip may cause inaccurate readings. It is important to not bend the test strip when attempting to apply the sample.

Tip Three - Using a Sideways Approach

When testing, it may be helpful to hold the meter on one side rather than flat or face up. In this position, the test strip sticks out of one end and the side with the buttons faces towards you.

The idea here is that the finger can be brought to the test strip from below until it touches the tip of the test strip. Then, the meter can be repositioned so the sample touches the absorbent opening at the end of the test strip. Notice how the meter can be positioned so that the test strip can point left or right. This allows for either hand to be used for testing.

Example

- **Step 1** - Gently massage the finger to bring blood to the puncture site.
- **Step 2** - Puncture the finger.
- **Step 3** - Hold the meter in your hand. Place your index finger under the test strip to help you guide the meter to where you have lanced.

- **Step 4** - Using overlapping dragging motions in an upward direction, repeat the motion around the puncture sight until you hear the meter beep then say, "testing."
- **Step 5** - Hold the test strip in place until you hear the meter say your test results.

Tip Four - Practice, Practice, Practice

Practicing is the key to feeling comfortable and confident in performing any new method, and we suggest using Prodigy Control Solution when practicing.

Practice until you can easily get a test result that falls within the target range for the Prodigy Control Solution three (3) times in a row. It may be helpful to have someone with sight watch you practice so they can provide feedback or answer questions you may have.

Applying Blood to the Test Strip

STEP 1: Insert the Test Strip Insert a test strip with the contact bar entering the test strip slot first. Push the test strip in as far as it will go without bending it. The meter will turn on automatically. See Figure 31.

Voice Speaks: "Your meter is on. Please wait. Please apply blood into the test strip."


Figure 31

STEP 2: Apply Blood Sample

When the meter shows "👉", apply blood to the open absorbent channel of the test strip where it will be drawn into the test strip. See Figure 32.

If the test strip confirmation window is full, you will hear a beep.

Voice Speaks: "Testing."


Figure 32

If the confirmation window is not full, “L-b” will be displayed. You will have to try the test again using a new test strip. See Figure 33.

Figure 33

If you have trouble filling the test strip, contact Customer Care at 1.800.243.2636 for assistance.

CAUTION

DO NOT:

- Smear or scrape the blood onto the test strip.
- Apply blood to the test strip when the test strip is not in the meter.
- Put blood or foreign objects into the test strip slot.

STEP 3: Read Your Result

After the meter counts down from 6 to 1, your blood glucose test result will appear, along with the unit of measure, date and time. See Figure 34.

Figure 34

Voice Speaks: “Blood glucose, 108 mg/dL.”

PRODIGY AutoCode® Blood Glucose Meter

Repeat Feature: After the test result is spoken, you may press the “**M**” button to hear the result again.

The blood glucose test result is automatically stored in the meter’s memory. Turn the meter off by removing the test strip. Discard the used test strip appropriately.

Important: If you do not apply a blood sample within three (3) minutes, the meter will automatically turn off. You must remove the test strip and re-insert it in order to turn on the meter and restart the testing procedure.

Caution: If you cannot test due to a problem with your testing supplies, contact Customer Care at 1.800.243.2636. Failure to test could delay treatment and/or lead to serious medical conditions.

Using the Meter Memory

The Prodigy AutoCode® meter stores a maximum of 450 of the most recent blood glucose test results with the date and time. It also provides 7, 14 and 28-day blood glucose test averages. You can review the individual or average test results by entering the memory mode.

STEP 1: Enter the Memory Mode

While the meter is turned off, press and hold the “**M**” button to turn the meter on. Press the “**M**” button again to enter memory mode. The 7-day average will appear, indicating that you are in the memory mode.

If you continue to press the “**M**” button, the 14 and 28-day averages will appear in order. You can then review the last 450 individual test results in the memory.

When using the meter for the first time or when the memory has been deleted, “- - -” appears, indicating there are no test results in the meter memory. See Figure 35.


Figure 35

STEP 2: Recalling Average Test Results

The 7-day average is calculated from the blood glucose results obtained during the last 7 days. It also indicates how many blood glucose tests have been performed within this period; e.g., 21 (21 tests in the last 7 days). See Figure 36.


Figure 36

The 14-day average is calculated from the blood glucose test results obtained during the last 14 days. See Figure 37.


Figure 37

It, too, indicates how many blood glucose tests have been performed; e.g., 41 (41 tests in the last 14 days). The 28-day average shows the same information.

STEP 3: Recalling Individual Test Results

After the 28-day average, the most recent test result with date and time will be shown. Press the "M" button once and the next most recent test result will appear. See Figure 38.

Each time you press and release the “M” button, the meter will recall up to your last 450 test results in order.

When the memory is full, the oldest test result will be deleted and the newest will be added. After reaching the last individual result, press the “M” button and the meter will turn off.


Figure 38

STEP 4: Exit the Memory Mode

At any point, you may press and hold the “M” button for four (4) seconds to exit memory mode.

Voice Speaks: “Please insert the test strip.”

Important: If no buttons are pressed for one (1) minute, the meter will display “OFF” and turn off automatically.

Viewing Results on a Personal Computer

Test results in your meter’s memory can be transmitted to your personal computer. Please visit the Prodigy web-site, www.prodigymeter.com, for more information.

SPECIFICATIONS

Dimensions: 3.77 in. (L) x 2.11 in. (W) x 0.86 in. (H)
95.7 mm (L) x 53.7 mm (W) x 21.8 mm (H)

Weight: 2.78 oz with batteries (79g)

Power Source: Two (2) 1.5 Volt AAA Alkaline Batteries

Battery Life: Over 1,000 tests

Display: LCD

Memory: 450 test results with date and time

External Output: USB Data Port

Auto electrode inserting detection

Auto sample loading detection

Auto reaction time count-down

Temperature warning

Operating Condition: 50°F-104°F (10°C-40°C),
10-85% R.H. (non-condensing)

Storage/Transportation Condition:

39.2°F -104°F (4°C-40°C), 10-85% R.H.

Measurement Units: mg/dL or mmol/L

Measurement Range: 20- 600 mg/dL (1.1-33.3 mmol/L)

The specifications may change without prior notice.

Caring for your Prodigy AutoCode® Meter

Follow the guidelines below to properly clean and disinfect your system:

- Please clean and disinfect your Prodigy AutoCode® meter at least once per week.
- When disinfecting your meter, use the PDI Super Sani-Cloth Germicidal Disposable Wipe or other disinfectant wipes with the same EPA number (#9480-4).
- If the wipe is dripping, gently squeeze out excess liquid before wiping the meter.
- Never immerse the meter in any type of liquid.
- The total time needed for the disinfectant solution to remain wet on the surface of the meter is two (2) minutes to ensure that total disinfection of all stated label claims are achieved.
- When cleaning and disinfecting the meter, be extremely careful to avoid letting any liquid enter the meter's test strip port or your mobile device's data or headphone ports, as this can cause inaccurate test results.
- We do not advise wrapping the meter in a wipe. Instead, wipe ALL exterior surfaces of the meter per the product label instructions. Be sure that the minimum contact time is met.

- Let meter air dry before testing or storing.
- Dispose of used wipes after cleaning and disinfecting.
- For more information on the PDI Super Sani-Cloth Wipes, visit www.pdipdi.com.
- PDI Super Sani-Cloth Wipes can be purchased at Walmart (www.walmart.com), Amazon (www.amazon.com) and many other online retailers, e.g. Moore Medical (<http://www.mooremedical.com>), Bio-medical (<http://bio-medical.com>), Quick Medical (<http://www.quickmedical.com>), etc.
- Always wash hands thoroughly with soap and water after handling the meter, lancing device or test strips.
- If meter deterioration occurs (for example, foggy screen, failure to detect test strip, faded or missing data on screen, false high or low readings), stop using the meter and contact Customer Care at 1.800.243.2636.

Cleaning your System

Cleaning is a process of removing dust, dirt, blood or other visible debris from the meter.

You should clean your Prodigy AutoCode® meter whenever it is visibly dirty by wiping the outside of the meter using a PDI Super Super Sani-Cloth Germicidal Disposable Wipe. Do not

use bleach, harsh abrasives, or cleaning materials other than the PDI Sani Cloths. The meter has been tested to withstand cleaning for up to two times per week for 3 years (520 validated cycles (1 cycle = 1 cleaning wipe and 1 disinfection wipe)).

Do not get any liquids, dirt, dust, blood, control solution, or cleaning materials inside the meter through the test strip port or the data port. Never spray cleaning solution on the meter or immerse it in any liquid.

Lancing Device

Wipe your Prodigy lancing device with a PDI Super Sani-Cloth Germicidal Disposable Wipe. Do not immerse the lancing device in any liquid. Do not use bleach. Disinfect your Prodigy lancing device twice a week.

After cleaning the device, unscrew the cap; wipe it with a PDI Super Sani-Cloth Germicidal Disposable Wipe for two (2) minutes. Make sure the cap is completely dry before reattaching it to the device.

Disinfecting your System

Disinfection is a process that destroys pathogens, such as viruses and other microorganisms on the meter and lancing device surfaces.

To disinfect:

1. The meter should always be cleaned before you disinfect it.
2. Take a new PDI Super Sani-Cloth Germicidal Disposable Wipe and wipe the entire surface of the meter and keep it wet for two minutes using the PDI wipe. After two minutes let the meter air dry.

If your meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be cleaned and disinfected prior to use by the second person. When the meter is only handled by one person, cleaning and disinfection should be performed periodically, at least once a week.

Storing Your System

Store the meter, test strips and control solution in the carrying case after each use in a cool, dry place.

Do not freeze.

Do not store in the kitchen or bathroom.

*Your meter is a precision instrument. Please handle it with care.

Battery

Your Prodigy AutoCode® meter comes with two (2) 1.5V AAA alkaline batteries. Batteries are manufactured in a charged state and not designed for recharging. Recharging the batteries can cause battery leakage, or in some cases, high pressure rupture. **Replace both batteries at the same time.**

 **Warning:** Batteries can explode or leak and cause burns if installed backwards, disassembled, charged or exposed to water, fire or high temperatures.

Low Battery

The meter will alert you when the power is getting low by displaying two (2) different messages:

1. When the " " symbol appears alone on the display, the meter is functional and the results remain accurate, but you should change the batteries as soon as possible.
2. When the " " symbol appears together with the "E-b" symbol on the display, the battery does not have enough power for a test. You must change the battery before using the meter.

Battery Replacement

When replacing the batteries, make sure the meter is turned off.

STEP 1: With the meter off, press the buckle on the battery cover and lift up to open the cover.

STEP 2: Remove the old batteries and replace with new ones. Be sure to align the batteries properly.

STEP 3: Close the battery cover. If the meter does not power on after you have replaced the batteries, check that the batteries are correctly installed.

STEP 4: Turn the meter on by pressing the **"M"** button to check the time and date. Replacing the batteries does not affect the meter's memory (previous test results stored in the memory). However, the date, time and unit of measurement settings may need to be updated.

Caution: As with all small objects, the batteries should be kept away from small children. If the battery is swallowed, seek medical assistance immediately.

Remove the batteries if you are not going to use the device for a long period of time (i.e., 3 months or more), or the batteries might leak chemicals.

System Troubleshooting

Special Messages

Special messages and error messages help to identify certain problems but do not appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without returning an error message.

In the event of a problem, refer to the information under “Action” in the “Error Messages” section. If you continue to have a problem, please refer to the “Troubleshooting Guide” section. If you follow the actions recommended but the problem is not resolved, please contact Customer Care at 1.800.243.2636 for assistance.

Message

What it Means

“Lo” appears when your result is below the measurement limit, which is less than 20 mg/dL (1.1 mmol/L).

“Lo” indicates hypoglycemia (low blood glucose). You should immediately consult with your healthcare professional.

Message	What it Means
	<p>"HI" appears when your result is above the measurement limit, which is higher than 600 mg/dL (33.3 mmol/L). You should immediately consult with your healthcare professional.</p>
Error Message	Description
	<p>What it means: Appears when the batteries cannot provide enough power for a test. Action: Replace the batteries immediately.</p>
	<p>What it means: The strip is used or moistened. Action: Test with a new test strip. If the problem persists, please contact Customer Care at 1.800.243.2636.</p>

Error Message	Description
	<p>What it means: Appears when the temperature is out of the system operating range (50°F-104°F or 10°C-40°C).</p> <p>Action: Repeat the test after the meter and test strip are within the operating temperature range.</p>
	<p>What it means: Insufficient blood sample. Test strip may be absorbing blood too slowly.</p> <p>Action: Insert a new test strip. If the problem persists, please contact Customer Care at 1.800.243.2636.</p>

All error messages are listed.

Troubleshooting Guide

- ❖ **The meter does not display a message after inserting a test strip**

Probable Cause	Actions
Battery exhausted.	Replace the batteries.
Battery incorrectly installed or absent.	Check that the batteries are present and correctly installed.
Test strip inserted upside down or incompletely.	Insert the test strip correctly with the contact bar entering the test strip slot.
Defective meter.	Please contact Customer Care at 1.800.243.2636 for assistance.

- ❖ **The test does not start after applying the sample.**

Probable Cause	Actions
Insufficient blood sample.	Repeat the test using a new test strip and a larger blood sample.

- ❖ The test does not start after applying the sample.

Probable Cause	Actions
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic shut-off. (Three (3) minutes after user action)	Repeat the test with a new test strip. Apply sample only when the "  " symbol appears on the display.
Defective meter.	Please contact Customer Care at 1.800.243.2636 for assistance.

❖ **The control solution test is out of range.**

Probable Cause	Actions
Error in performing the test.	Read the instructions thoroughly and repeat the test again.
Control solution bottle not shaken well.	Shake the control solution bottle vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiration date and the discard date of the control solution.
Control solution that is too warm or too cold.	Control solution should come to room temperature before testing.
Test strip deterioration.	Repeat the test with a new test strip.

❖ The control solution test is out of range.

Probable Cause	Actions
Meter malfunction.	Please contact Customer Care at 1.800.243.2636 for assistance.

Information About Your Prodigy AutoCode® Meter

Comparing Meter and Laboratory Results

The test results you obtain from your meter may be different from your laboratory results due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way (see test strip package insert for typical accuracy and precision data as well as important information on limitations). To make an accurate comparison between meter and laboratory results, follow the guidelines below.

Before you go to the laboratory:

- Perform a control solution test to make sure that the meter is working properly.
- It is strongly recommended to fast at least eight (8) hours before doing comparison tests.
- Take your meter with you to the laboratory.

While at the laboratory:

- Make sure that the samples for both tests (the meter test and the laboratory test) are taken and tested within 15 minutes of each other.
- Wash your hands before obtaining a blood sample.

- Never use your meter with blood that has been collected in a gray top test tube.
- Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods, especially if you have recently eaten, exercised, taken medication or experienced stress.¹ In addition, if you have eaten recently, the blood glucose level from a finger stick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test.²

Therefore, it is best to fast eight (8) hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of bodily fluids (severe dehydration) may also cause a meter result to be different from a laboratory result.

References: ¹Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51. ²Sacks, D.B.: "Carbohydrates." Burtis, C.A, and Ashwood, E.R. (ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

Performance Characteristics

- **Accuracy:** > 95% of the individual glucose results fell within ± 15 mg/dL (0.83 mmol/L) at glucose concentration < 75 mg/dL (4.2 mmol/L) and within $\pm 20\%$ at glucose concentration ≥ 75 mg/dL (4.2 mmol/L), compared with glucose analyzer YSI 2300.

Table 1 represents samples for glucose results lower than (<) 75 mg/dL

Difference between laboratory and Prodigy AutoCode® meter	Within ± 5 mg/dL	Within ± 10 mg/dL	Within ± 15 mg/dL
number (%) for Finger Blood	9/14 (64%)	13/14 (93%)	14/14 (100%)
number (%) for palm blood	5/14 (36%)	14/14 (100%)	14/14 (100%)
number (%) for forearm blood	3/14 (21%)	13/14 (93%)	14/14 (100%)
number (%) for upper arm	8/14 (57%)	14/14 (100%)	14/14 (100%)
number (%) for calf blood	6/14 (43%)	12/14 (86%)	14/14 (100%)
number (%) for Thigh blood	4/14 (29%)	13/14 (93%)	14/14 (100%)

PRODIGY AutoCode® Blood Glucose Meter

Table 2 represents samples for glucose results greater than (\geq) 75 mg/dL.

Difference between laboratory and Prodigy AutoCode® meter	Within $\pm 5\%$	Within $\pm 10\%$	Within $\pm 15\%$	Within $\pm 20\%$
number (%) for Finger Blood	30/86 (35%)	59/86 (69%)	85/86 (99%)	86/86 (100%)
number (%) for palm blood	29/86 (34%)	54/86 (63%)	83/86 (97%)	86/86 (100%)
number (%) for fore-arm blood	28/86 (33%)	59/86 (69%)	83/86 (97%)	86/86 (100%)
number (%) for upper arm	32/86 (37%)	63/86 (73%)	82/86 (95%)	86/86 (100%)
number (%) for calf blood	30/86 (35%)	53/86 (62%)	84/86 (98%)	86/86 (100%)
number (%) for Thigh blood	23/86 (27%)	48/86 (56%)	83/86 (97%)	86/86 (100%)

Linear regression analysis:

Finger: $y = 1.0177x - 2.6955$; $R^2 = 0.975$

Palm: $y = 1.0437x - 6.3933$; $R^2 = 0.9698$

Forearm: $y = 0.9628x + 5.9635$; $R^2 = 0.9675$

Upper arm: $y = 1.0175x - 3.9439$; $R^2 = 0.9681$

Calf: $y = 1.0058x - 2.4028$; $R^2 = 0.9684$

Thigh: $y = 1.0137x - 3.3484$; $R^2 = 0.9576$

- **Precision:** This study shows the CV (correlation variation) is less than 5%.

Expected Test Results

Please consult your healthcare professionals to establish your own target ranges. The normal adult fasting blood glucose range for a nondiabetic person is less than 100 mg/dL (5.55 mmol/L) and less than 140 mg/dL (7.77 mmol/L) up to 2 hours after meals.

*Source: American Diabetes Association. Classification and diagnosis of diabetes. Sec. 2. In Standards of Medical Care in Diabetes - 2016. Diabetes Care 2016;39(Suppl. 1):S16.

The Prodigy AutoCode® meter is designed to help you and your healthcare professional manage your diabetes. You must always rely on your healthcare professional to interpret your test results and to decide how to treat your diabetes.

Unexpected Test Results

Unexpected test results can occur. When this happens, please refer to these cautions.

Low Blood Glucose Test Results

False low results may occur if you are severely dehydrated. If you think you are severely dehydrated, contact your healthcare professional immediately.

If your test result is lower than 70 mg/dL or “Lo” appears on your meter screen, this could be a result of hypoglycemia. Although this result could be due to a testing error, it is best to treat first and retest later. Contact your healthcare professional immediately. See Figure 39.

Figure 39

High Blood Glucose Test Results

If your test result is higher than 180 mg/dL, this could be a result of hyperglycemia. If you are uncertain about this test result, retest. If the result continues to be higher than 180 mg/dL, contact your healthcare professional immediately.

If “HI” appears on your meter screen, this could be a result of severe hyperglycemia (results higher than 600 mg/dL). Retest your blood glucose. If “HI” appears again, contact your healthcare professional immediately. See Figure 40.

Figure 40

Repeated Unexpected Test Results

If you continue to get unexpected test results, check your blood glucose monitoring system by performing a control solution test. Refer to “Performing a Control Solution Test” section in this manual.

Unusual Red Blood Cell Count

A hematocrit range (percentage of red blood cells in the blood) that is extremely high (above 60%) or low (below 20%) can also cause false readings.

Accuracy of Meter Measurements

FAQ:

Can two (2) blood glucose meters produce different readings that are from the same blood sample? Does this make one reading more accurate than the other?

Answer:

Yes, different meters can produce different readings using the same blood sample.

No, this does not necessarily mean that one result is more accurate than the other.

The reasons for this are as follows:

1. The FDA recognized standard for glucose meter accuracy is ISO 15197: In Vitro Diagnostic Test Systems Requirements for Blood Glucose Monitoring Systems for Self Testing in Managing Diabetes Mellitus. This standard sets accuracy requirements that manufacturers must meet before being approved by the FDA to market their glucose meters.

2. The accuracy requirements for the blood glucose monitoring system—meters and test strips—as stated in ISO 15197 are:

Ninety-five percent (95%) of the individual glucose results shall fall within ± 15 mg/dL (0.83 mmol/L) of the results of the manufacturer's measurement procedure at glucose concentrations < 75 mg/dL (< 4.2 mmol/L) and within $\pm 20\%$ at glucose concentrations ≥ 75 mg/dL (≥ 4.2 mmol/L).

How does this answer the FAQ?

The accuracy of the system is determined by a clinical study using 100 blood samples at interval glucose concentrations ranging from < 50 mg/dL to over 400 mg/dL measured using both the manufacturer's meter and a clinical blood glucose analyzer for comparison (reference measurement).

Example

If the manufacturer's reference measurement is 70 mg/dL this means that the meter measurements are considered accurate if they fall within ± 15 mg/dL of 70 mg/dL (from 55 mg/dL to 85 mg/dL).

If meter A produces a measurement of 60 mg/dL and meter B measures 80 mg/dL, then both are in the acceptable range and meet the accuracy requirement.

Each manufacturer has met the System Accuracy requirement as part of FDA testing to be cleared for marketing. In the System Accuracy testing, a series of blood glucose readings will result in various glucose concentration levels for one manufacturer's system. This is the same for other manufacturers' systems. Because of the range of acceptable readings for the accuracy, identical readings from two different meters on the same blood sample may not happen.

Symbol Information

Symbol	Referent
	Do not re-use. Single use only.
	Consult Operating Instructions.
	Keep away from sunlight.
	Keep dry.
	Temperature limitation.
	Use by.
	Date of manufacture.
LOT	Batch code.
REF	Catalog number.
SN	Serial number.
CONTROL	Control.

Warranty Information

Limited Lifetime Warranty

Prodigy Diabetes Care, LLC extends a limited lifetime warranty to consumers who purchase a new Prodigy brand meter. Under this limited lifetime warranty, your new meter is covered from the original date of purchase as long as it has not been modified, altered, or misused, and Prodigy will replace the meter free of charge if it is defective in material or workmanship.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE. PRODIGY DIABETES CARE, LLC WILL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING DIRECTLY OR INDIRECTLY FROM THE FAILURE OF THE PRODUCT TO PERFORM IN ACCORDANCE WITH SPECIFICATIONS.

Some states do not allow the exclusion or limitation of other express or implied warranties or incidental or consequential damages, so the above limitations or exclusions may not apply to you.

PRODIGY
diabetes care

Prodigy Diabetes Care, LLC

PO Box 481928

Charlotte, NC 28269

1.800.366.5901

Customer Care:

1.800.243.2636

Monday - Friday 8:30AM to 5:30PM (EST)


www.prodigymeter.com

©2017 Prodigy Diabetes Care LLC. Prodigy®, AutoCode®, and Prodigy AutoCode® are registered trademarks of Prodigy Diabetes Care, LLC. All other trademarks are property of their respective owners. All rights reserved. 18-047 ACOM01 Rev.1 8/18

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