
# Molift EvoSling LowBack


User manual

**molift**<sup>®</sup>  
a part of Etac

BM28499 Rev. A 2017-01-26

CE

## Contents

User manual - English - 4
Bruksanvisning - Svenska - 8
Brukermanual - Norsk - 12
Brugsvejledning - Dansk - 16
Käyttöohje - Suomi - 20
Gebrauchsanweisung - Deutsch - 24
Handleiding - Nederlands - 28
Manual d'utilisation - Français - 32
Manual utente - Italiano - 36
Manual de usuario - Español - 40

**General**

**Declaration of conformity**  
 The Molift EvoSling LowBack and related accessories described in this operator manual are CE marked in accordance with EU Council Directive 93/42/EEC concerning medical devices, class 1, and has been tested and approved by a third party according to standard EN ISO 10535:2006.

**Conditions for Use:** Lift and transfer of a person will always pose a certain risk and only trained personnel should use the equipment and accessories covered by this user manual.

**Warranty:** 2-year warranty against defects in workmanship and materials of our products. For Terms and conditions, see [www.etac.com](http://www.etac.com)  
 Visit [www.etac.com](http://www.etac.com) to download the current documentation.

**Important**

This User Manual contains important safety instructions and information regarding the use of the sling and accessories. In this manual the user is the person being hoisted. The assistant is the person operating the hoist.

 This symbol indicates important information related to safety. Follow these instructions carefully.

 **Read User Manual before use!**  
 It is important to fully understand the content of the user manual before attempting to use the equipment. Read the manual for both hoist and sling

**General Safety Precautions:** Molift EvoSling LowBack is designed for use with a 2-point suspension but also works with a 4-point suspension. Only use accessories and slings that are adjusted to fit the user, type of disability, size, weight and type of transfer. It is important that the sling has been tested with the individual user and for the intended hoisting situation. Make a decision on whether one or more assistants are required. Plan the hoisting operation in advance to ensure that it is as safe and smooth as possible. Remember to work ergonomically.


**Assess the risk and take notes.**  
 You as the assistant are responsible for the safety of the user!

 Molift slings should only be used to lift people. Never use the sling to hoist or move objects of any kind.

 If maximum load (SWL) differs between the hoist, suspension and body support unit, then the lowest maximum load should always be used

 The sling should not be stored in direct sunlight.

**Labels and symbols**


1. Main label
2. Size label
3. Product label with barcode that contains:  
 (01)xxxxxxxxxxxxxxxx(21) xxxxxxxx  
 (01:EAN number (21):Serial number
4. Periodic inspection label / Name tag

- Symbols:**
-  - Date of manufacture. YYYY-MM (year/month)
  -  - CE marked
  -  - Refer to user manual
  -  - This way up, this side is outside
  -  - Max user weight


- Wash symbols:**
-  - Slings can be washed at temperatures ranging from 60 up to 85° C.
  -  - Tumble drying, max 60° C.
  -  - Do not iron
  -  - Do not dry clean
  -  - Do not bleach

**About EvoSling LowBack**

Etac supplies a wide selection of slings for different types of transfers. Molift EvoSling LowBack has a low back height, split leg support with a medium opening. The sling can be used for lifting to and from a sitting position. Approved for users weighing up to 300kg. Molift EvoSling LowBack is available in sizes XS – XXL, in padded polyester and net polyester with or without padded leg supports. The net version sling has a centre mark line and both versions have application pockets to facilitate application of the sling. The sling has a guide handle for positioning of the user when lifted. The user can have their arms inside or outside the sling. Molift EvoSling is developed to be combined with a 2-point suspension but the sling also works with a 4-point suspension. See the combination list for the correct sling and suspension.

**Technical data**

**Safe Working Load (SWL):**  
XXS-XL: 300 kg (660 lbs)

**Material:**  
Polyester, Polyethylene

Size	A	B
XS	41	40
S	43	42
M	44	48
L	47	55
XL	51	62
XXL	55	70


**Expected Lifetime:**

The product has an expected service life of 1 to 5 years under normal use. The service life of the product varies depending on usage frequency, materials, loads and how often it is washed. Washing at higher temperatures wears out the material faster.

**Lifting Strap Ladder loops**

The ladder loops on the lifting band give different options for mounting the sling to the suspension. Use the ladder that makes the sling most comfortable for the user.

The ladder loops on the straps have different colours so the strap can easily be fastened with equal length on left and right side. Distance between each ladder loop (A) is 9cm.

**Combination list**

 2 Point Suspension	Sling size		
	XS-S	M-L	XL-XXL
Small 340mm	✓	✓	
Medium 440mm	✓	✓	✓
Large 540mm		✓	✓

 4 Point Suspension	Sling size		
	XS-S	M-L	XL-XXL
Small 340mm	✓	✓	
Medium 440mm	✓	✓	✓
Large 540mm		✓	✓

**Before use / Daily check**

 **Before using the sling for first time it should be marked with the date of first time use on periodic inspection label.**

Inspection to be performed daily or before use:

- Make sure the user has the required ability for the sling
- Make sure the sling is suitable for the hoist/ suspension to be used.
- Make sure the sling does not have visible damage or frays.


**Periodic inspection**

Periodic inspection of the non-rigid body support unit should be undertaken at the time intervals stated by the manufacturer, but at least every 6 months. More frequent inspections may be required where a non-rigid body support unit is used or cleaned more frequently than normal.

 **Never use a faulty or damaged sling as it can break and cause personal injury. Destroy and discard damaged and old slings**

 **Make sure not to damage or remove labels when cleaning.**

**Size guide**


Measurements in the table should be used as a reference. The correct size depends on the the users width and back height when they are in a sitting position. Measure the width and back height of the user in a sitting position. Choose the sling with the most correct measurements to make sure it fits.



Size	User weight - kg (lbs)	Width (W) - cm	Height (H) - cm
XS	17-25 (37-55)	31-38	60-70
S	25-50 (55-110)	35-42	70-80
M	45-95 (99-210)	37-47	80-95
L	90-165 (198-353)	40-51	95-100
XL	160-240 (353-529)	45-70	100-105
XXL	230-300 (507-660)	50-75	100-120

**Sling checkpoints**

1. The suspension is not too close to the users head
2. Leg straps are crossed
3. The leg support is placed under the users thigh without creases.



### Sitting position - Sling application

1. Ask the user to lean forward, support them if necessary. Slide the sling down to the level of the tailbone/coccyx.

2. Pull the leg supports forward and position them under the user's thighs.



3. Make sure the lifting straps are the same length on both sides and that the sling is placed under the users thighs without any creases. Cross the leg loops by pulling one of the loops through the other and connect all 4 lifting straps to the suspension. Make sure the matching lifting straps are fastened with equal length on the left and right side. Make sure the user is in a comfortable and secure position in the sling.

1. Check that the sling is correctly fitted around the user and that the strap loops are correctly fitted to the suspension hooks.
2. Start lifting until the lifting straps are stretched without lifting the user. Ensure that all four loops of the sling are securely fastened to avoid the user slipping or falling out of the sling.
3. Hoist the user, and perform the transfer.

Be careful during movement, the user may swing during turns, stops and starts. Be careful when manoeuvring close to furniture and similar to ensure the suspended user does not hit anything.

Avoid deep pile carpets, high thresholds, uneven surfaces or other obstacles that may block the castors. The hoist may become unstable if forced over such obstacles increasing the risk of tipping over.

 User must be facing the hoist when lifting or lowering from or into a chair.


	
---	--

 Never leave a user unattended in a hoisting situation.

 The hoist should not be used to lift or move users on sloping surfaces


### Transfer

When moving the user, stand to the side of them. Make sure that their arms and legs do not obstruct the seat, bed, etc.

 Never hoist the user higher than necessary to carry out a lift. Remember that wheels on a mobile lift must NOT be locked.

### Lowering into sitting


When transferring the user into a sitting position, the following techniques can be used for lowering the user as far back into the seat as possible (Use A or B):


- A Push gently on the users knees to guide the user in towards the backrest of the chair.
- B Tilt the seat backwards.

 Do not hoist with the guide handles! Excessive force on the handles can cause the sling to tear.

## Sitting position - Removing the sling

1. Lower the user into a sitting position. Release the lifting straps from the suspension. Remove the leg support by folding the leg support under the sling itself and the users thigh and pull out.

2. Stand next to the chair and ask the user to lean forward. Support the user with one hand if necessary. Pull the sling away with one hand, and ask the user to lean back in the chair.

Remove the sling gently to ensure the user does not fall forward.

## Recycling

Recycling instructions is available at [Etac.com](http://Etac.com)

## Troubleshooting

Symptom	Possible Cause/Action
The sling does not fit the user.	Wrong sling size. Try another size.
	The sling is too small when the leg supports do not reach halfway around the inner thigh and the suspension is too close to the user's face. Try another size.
	The sling is not correctly positioned. Pull the sling further down at the back when fitting it.
All the straps do not lift at the same time.	If the sling straps lift unevenly, it means that the sling can slide up at the back which can lead to the user sliding down into the opening of the sling. Make sure that the leg supports are positioned correctly on the thighs.
User is sliding through/out of the sling.	The sling is not connected correctly according to the users individual prescription.
	The user is wearing "slippery" clothes. Change or remove some of the clothing

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